



Bulk Hippex Powder is delivered in 3 packings.

- 360 gram in large pots, to prepare up to 4 liter
- 1800 gram in a small container, to prepare up to 20 liter
- 4500 gram in a large container, to prepare up to 50 liter

The image shows the 4 different packings: the small jars containing 22.5 gram to prepare 250 ml extender, a pot containing 360 gram to prepare up to 4 liters, and 2 containers that have 1800 and 4500 gram for up to 20 and 50 liters, respectively.



The use of bulk powder has 3 main advantages:

1. The price of Hippex Bulk Powder is lower than the price of Hippex in small jars.
2. Bulk powder offers the possibility for end users to weigh the exact amount of Hippex powder needed. Then you can prepare the right amount of extender just before semen collection.
3. Distributors may make their own packings. This offers the possibility to pack in blisters, as frozen liquid medium, or as ready-made extender.

When you have warm water (distilled or demineralized) available, you can prepare the needed amount of extender in an instant. Just weigh the appropriate amount by calculation of the amount of powder to be weighed based on 90 gram Hippex powder per liter water. You may keep warm water in an incubator or water bath at all times, to have it available when needed.

It is recommended to store Bulk Hippex Powder in the refrigerator. When taking the powder from the refrigerator and thereafter weighing the appropriate amount of powder, work quickly to prevent moisture from the air precipitating on the cold powder.

This depends of course on the outside temperature and relative humidity. When too much moisture precipitates onto the powder, it will have in time a negative effect on the quality of the powder.

In a warm and humid environment (> 20 °C and > 50 % RH), it is best to take the Bulk Hippex Powder out of the refrigerator the day before preparing the extender. In this way, too much moisture precipitation onto the powder will be avoided. The powder will then also dissolve easier in water.

Each customer who orders at least the amount of powder to prepare 24 liter Hippex medium will receive also 2 spoons and a precision scale with 2 decimals (max. 350.00 gram). A scale is also a useful tool to measure (in grams) the amount of water needed.